IKATRON® ETS-D4 fuzzy

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OPERATING INSTRUCTIONS EN 9

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DE

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CE-DECLARATION OF CONFORMITY

EN

We declare under our sole responsibility that this product corresponds to the regulations 89/336/EG and 73/23/EG and conforms with the standards or standardized documents DIN EN IEC 61010-1 and DIN EN IEC 61326-1.

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In order to be able to use the appliance properly and safely, every user must first read the operating instructions and observe the safety instructions contained therein. Take care of these operating instructions and keep them in a place where they can be accessed by everyone.

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Warranty

You have purchased an original IKA laboratory machine which meets the highest engineering and quality standards.

In accordance with IKA guarantee conditions, the guarantee period is 24 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

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Safety instuctions

Note: With this equipment only media may be regulated, whose flash point lies over the adjusted safety temperature delimitation of the used magnet agitator and/or heating system.

The instrument may not be operated in rooms with explosion hazard.

The stainless steel temperature sensor may not be operated because of corrosion danger in aggressive media such as acids or caustic solutions or distilled water. Use for this the glass sensor.

Please make shure that the spiral cable does not contact the heating plate.

Caution when affecting the temperature sensor!

Burn danger! The temperature sensor can become hot in the medium up to 400°C !

Observe the accident prevention requirements applicable to the activity and the work station.

Wear personal safety equipment.


Danger:

In the case of use the IKATRON® ETS-D4 fuzzy with heaters not over a connection after DIN 12878 and thus no break fuse have, must be examined before start up the normal and safe enterprise by the user.

(A break fuse is a mechanism guarantees the fact that with break of the contact thermometer (ETS-D4) and with the separation of the connection of the heating electric circuit is switched off).

Corect use

The IKATRON ETS-D4 fuzzy is an instrument for the precise control of temperatures. The sensor of the instrument must be immersed at least 20mm deep in the medium to be tempered. The instrument can be connected to any magnetic agitator with contact thermometer connection which conforms with the requirements listed in the Technical Data (see section Technical Data).

Unpacking

Please unpack the machine carefully and inspect for damage. It is important that any transit damage should be noted at the time of unpacking. In certain circumstances it may be necessary to investigate immediately (post, rail or freight forwarder).

As supplied, the package contains: one IKATRON ETS-D4 fuzzy and one set of operating instructions.

Factory settings

The ETS-D fuzzy with the following factory settings one delivers:

Mode of operating	A
Set-point temperature	-10 °C
HI-TEMP	400 °C
Mode of thecontroller	Fuzzy-Controller
Type of the sensor	H 62
Extension cable	---
Error 5 response time	3 min

Put back the parameter to factory settings (RESET)

For putting the parameters back to factory settings the following steps are necessary:

- Switch off the equipment
- The foil keys “+” and “-” keep pressed and the equipment switch on after approx. 5 secondes flashes the announcement, now the foil keys release
- Now the equipment is again neglected to factory settings (see above)

Commissioning

Plug the instrument into the contact thermometer jack on the magnetic agitator (e.g. IKA RCT basic or IKA RET basic). The instrument is activated when the agitator is switched on. If there is no voltage supply between 8V and 16V DC at the contact thermometer jack, the instrument must be run off the H52 power supply kit (H50 twin-jack adapter and H51 plug power supply). The instrument will then receive the required voltage from the power supply or twin-jack adapter.

An automatic self-test is performed when the instrument is activated. During the self-test all the segments of the LCD display are visible (**Fig.1**). Then the LCD displays illustrated in **Fig.2** to **Fig.7** appear (standard factory setting).

System configuration

Before starting to work with the ETS-D4 fuzzy the system configuration must be set. To do this, hold down the “-” key when switching on.

1. Setting the type of sensor used

The LCD display indicates the set type of sensor, e.g. H 62 (**Fig.8**). Accessories see page 7. Use the “+” and “-” keys to set the type of sensor H 60 ... H 67 you are using. To change the setting, hold down one of the two keys for approx. 2 seconds. The displayed sensor is set and stored if no key is pressed for approx. 5 seconds.

2. Setting the extension cable used

The LCD display indicates the set type of extension cable (**Fig.9**) in small letters. The setting is made by pressing the “+” and “-” keys. Possible settings are: H 70, H 71 or — if no extension cable is used. The displayed extension cable is set and stored if no key is pressed for approx. 5 seconds.

3. Setting the control mode

The LCD display indicates the control setting. Three control modes are available. The setting can be changed by pressing the “+” and “-” keys. The displayed control mode is set and stored if no key is pressed for approx. 5 seconds.

- a) Fuzzy control (**Fig.10**). Standard factory setting, best control behavior.
- b) Fast Fuzzy control (**Fig.11**). Faster warming up than with a), very good control behavior in stabilized state. “>” displays Fast Fuzzy control mode.
- c) Two-Point control (**Fig.12**). Faster warming up than with a) and b), severe overswinging on first warm-up.
Possible application: Using the ETS-D4 fuzzy with KSR 11 as a limit monitor. The heater or instrument remains switched on until the set nominal temperature is exceeded. “:” signifies Two-Point control mode.

4. Adjustment of the Er 5 response time

This is the amount of time required to detect that the measurement sensor is not in the medium (when the heater is turned on).

Safety information: This detection is not activated unless the difference between the sensor temperature and the target temperature is >5K if the measurement sensor temperature remains constant (+/-0,5K).

The adjusted Er5 response time is displayed in the LCD display, for example 3minutes (**Fig.4a**). You can change the setting to any value between 3 and 20 minutes (in increments of minutes) by pressing the "+" or "-" key. The adjusted time is accepted by the system if no key is pressed for about 5 seconds.

Operating modes

The ETS-D4 fuzzy can be run in three different operating modes. (Displays see **Fig.13** to **Fig.15**)

Operating mode "A" The instrument is factory-set to this mode. When the instrument is switched on the nominal temperature is always set to -10°C. The nominal temperature can be adjusted with the "+" and "-" keys. The HI-TEMP (maximum nominal temperature setting) can only be set in this mode.

Operating mode "B" In this mode the most recently set nominal temperature is adopted and becomes active when the instrument is switched on again. The setting can be adjusted with the "+" and "-" keys.

Operating mode "C" In this mode the current settings are adopted and become active when the instrument is switched on again. The setting cannot be adjusted. Holding down the "+" key when switching on the instrument switches to the next operating mode in the sequence A-B-C-A-B-C-A etc.

Setting the HI-TEMP

The "HI-TEMP" (maximum nominal temperature setting) can only be set in operating mode "A". After the self-test, HI-TEMP appears in the LCD display (**Fig.16**). The HI-TEMP can now be set between 0 and 400°C using the "+" and "-" keys. The displayed HI-TEMP is set and stored if no key is pressed for approx. 5 seconds. The LCD display then shows "HI OK" (**Fig.17**).

Setting the nominal temperature

The nominal temperature is adjusted in 1K steps by pressing the "+" and "-" keys. If the key is held down for longer than 5 seconds, the nominal temperature setting changes in 10K steps. The LCD display shows: the current actual temperature with physical units, the nominal temperature at top right, and OFF to indicate that the heater of the connected device is switched off. (**Fig.18**)

If the set nominal temperature is higher than the current actual temperature, the word OFF in the display changes to ON and the heater of the connected device is switched on. (**Fig.19**)

The star (Q) in the LCD display indicates that the instrument is being operated in mode "B" or "C". (**Fig.20**)

Replacing extension cable and sensor

In order to replace the sensor or the extension cable, pull the sensor and the cap downward until the plastic snap fastener is released. Then push the sensor or extension cable with the cap over the fitting on the instrument to reestablish the connection.

Error messages

The following error messages may occur (**Fig.21**) and have the causes listed below:

- Er 1 Sensor not connected or defective (open circuit).
- Er 2 Analog output loaded with <math><50\text{k}\Omega</math>, output voltage too low.
- Er 3 Internal temperature above +76 °C.
- Er 5 Sensor not in medium (no change in medium temperature after heater has been on for 3 minutes).
- Er 6 Sensor short circuit.
- Er 7 The maximum set HI-TEMP is exceeded.
- Er 8 Calibration error (re-calibration by field service required).
- Er 9 IC107 (4066) defective
- Er 10 T_105 defective, OUT_TEMP >480 °C (4.8V)
- Er 11 Supply voltage <math><7.5\text{V}</math>
- Er 12 Supply voltage >16.1V

Analog output

The actual temperature is represented as a voltage (1V corresponds to 100°C) at the 3.5mm jack of the twin-jack adapter.

Tolerance: 0...100 °C ± 10 mV
 100...400 °C ± 25 mV

The internal resistance of the connected measuring instrument or recorder must be greater than 50 k Ω .

Accessories

H 36	Retaining bar		
H 44	Retaining clip		
H 50	Twin-jack adapter		
H 52	Power supply kit		
H 60	Short sensor	60 mm	stainless steel
H 62	Standard sensor	270 mm	stainless steel
H 64	Glass sensor	50 mm	glass
H 66	Glass sensor	270 mm	glass
H 70	Extension cable, 1 m long		
H 71	Extension cable, 2.5 m long		
Ak 2.5	Analog cable		

Technical data

Measuring range:	°C	-10 ... 400
Resolution:	K	0.1 (-10 ... +100 °C)
	K	1 (100 ... 400 °C)
Measuring accuracy:	K	± 0.3 + sensor tolerance Pt 1000 DIN IEC 751 Class B
Setting accuracy:	K	1
Control deviation:	K	± 1
Supply voltage:	V DC	8 ... 16
Current consumption:	mA	16 (at 9V)
Ambient temperature:	°C	0 ... 60
Relative humidity:	%	80
ON time:	%	100
Degree of protection:		IP 44
Cable length:	mm	350
Plug:		6-pin DIN 45 322
Connection:		DIN 12 878 Class 2
Dimensions: (W x D x H)	mm	70 x 150 x 65 (without sensor)
Weight:	kg	0.2

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
List of spare parts

Item	Designation
1	Rear wall
2	Clamp
3	Sealing ring
4	Spiral cable
5	Sealing insert
6	Union nut
7	Housing
8	Dished washer
9	Membrane keypad
10	Sensor jack; complete
12	BLP temperature regulator
13	Plastite screw 2.5 x 8

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IKA® Works, Inc.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

North America

IKA® Works, Inc.

2635 North Chase Pkwy SE
Wilmington, NC 28405-7419
USA
TEL. +1 800 733-3037
TEL. +1 910 452-7059
FAX +1 910 452-7693
E-mail: usa@ika.net

IKA® Works, (Asia) Sdn Bhd

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Asia - Australia

IKA® Works (Asia) Sdn Bhd

No. 17 & 19, Jalan PJU 3/50
Sunway Damansara Technology Park
47810 Petaling Jaya
Selangor, Malaysia
TEL. +60 3 7804-3322
FAX +60 3 7804-8940
E-mail: sales@ika.com.my

IKA® Japan Y.K.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Japan

IKA® Japan Y.K.

293-1 Kobayashi-cho
Yamato Koriyama Shi
639-1026 Japan
TEL. +81 74358-4611
FAX +81 74358-4612
E-mail: japan@ika.de

IKA® Works Guangzhou

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

China

IKA® Works Guangzhou

173-175 Friendship Road
Guangzhou Economic & Technological
Development Zone
Guangzhou 510730, P.R.CHINA
TEL. +86 20 8222-6772
FAX +86 20 8222-6776
E-mail: sales@ikagz.com.cn